# RaySafe X2

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## CONNECT TO COMPUTER USING BLUETOOTH

Connect the Bluetooth adapter to a SENSOR port on the base unit. A gray Bluetooth status symbol appears in the upper right corner of the screen.

X2 View will automatically search for base units. The first time you connect, you will have to select your base unit from the Bluetooth menu in X2 View. The Bluetooth status symbol will turn white when connected.

When restarting X2 View, the base unit will automatically get connected.

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## MEASURE WITH CT SENSOR

Push the connected sensor tight into the phantom ...

... or use the flexi stand for positioning free in air.

Expose.


Read the result.

**Note!** The effective length of the CT sensor is 100 mm, and is marked with center and edge lines on the tube.

**Tip!** Tap a parameter to get large digits, parameter information, and, if applicable, waveform. Switch between these by swiping sideways.

## MEASURE KVP ON CT

For kVp measurements, use an R/F sensor. Position it flat on the table and expose without rotation, using scout, tomogram or topogram mode. Use no or slowest table movement.

# CT SENSOR: MEASUREMENT PARAMETER DEFINITIONS

**Dose** is calculated from all recorded data.

**Dose rate** is the average dose rate, calculated as dose/time. For measurements longer than 3 s, intermediate readings are moving averages.

**Time** starts the first time the dose rate waveform reaches 50% of peak, and ends the last time it drops below 50%. Intermediate readings are time since start trig.

**Dose length product** is dose·10 cm.

## SETTING: DOSE UNIT

Select if all dose units should be displayed in Gray or Röntgen.


1 Gy = 114.1 R

X2 systems approved by PTB do not have the possibility to change the dose unit to Röntgen (R).


## SETTING: STOP DELAY

The stop delay setting defines for how long the instrument will wait for more radiation to include in the same measurement.

Use a longer stop delay when using a phantom and a slow rotation rate, to prevent that an axial scan will be cut into two measurements.


Stop delay ( $t$ ) setting longer than the time between pulses ( $t > T$ ), will result in one long measurement. Note that you have to wait as long as the stop delay setting ( $t$ ) after the last exposure before the measurement ends and values will be displayed.

Stop delay setting ( $t$ ) shorter than the time between pulses ( $t < T$ ), will result in many short measurements (one for each pulse).

## MEASURE WITH LIGHT SENSOR

**Luminance:** Place the connected sensor with the aperture centered on the area you want to measure, and with luminance selected on the sensor wheel.

**Illuminance:** Place the connected sensor with the diffuser in the direction you want to measure, and with illuminance selected on the sensor wheel.

The sensor is constantly measuring. If you want to store a value, press the button once the reading is stable. The measurement will continue below the stored value.

**Tip!** Tap a parameter to get large digits and parameter information. Switch between these by swiping sideways.

**Zero adjust:** If prompted to make a zero adjust, select zero adjust (0) on the sensor wheel. The zero adjust will take approximately 10 s.

**Note!** Store the sensor with the sensor wheel in a closed position (0 or illuminance), to prevent dust in the optics.

# LIGHT SENSOR: MEASUREMENT PARAMETER DEFINITIONS

**Illuminance** is the amount of light incident on a surface.



**Luminance** is the amount of light emitted from a surface.

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© 2014.03 Unfors RaySafe 5001088-3